

Section B **Audiometry Audiometers****Amplivox 160 Screening Audiometer**

The screening audiometer of choice for primary care and school use:

- Reduced test time with threshold retention function
- Compact and totally portable
- Supplied with mains adaptor, carrying case and patient response switch
- Archive quality print-out available with specified printer (optional)
- Easy and straightforward to use
- 10 comprehensive test frequencies
- -10 to 100 dBHL intensity range
- Includes internationally approved earphones
- UK design and manufacture
- Mains 240v powered



**AU160** Amplivox 160 Screening Audiometer

**£790.00**

**Optional accessories for Amplivox 160:**

<b>AU160.1</b> Battery power function (enables unit to be used with 4 x C cell batteries)	<b>£98.00</b>
<b>AU160.3</b> Printer and cable	<b>£245.00</b>

Audiometers / Tympanometers **Audiometry Section B****Interacoustics MT10 Tympanometer****Quick, lightweight tympanometry**

Built to accommodate modern healthcare settings where fast and reliable testing is paramount. Lightweight and incorporating an extended pressure range for differential diagnosis. Easy to use with large clear high resolution display.

- +300 to -600 daPa
- Charger base and integral printer
- RS232 interface
- 20 patient memory
- NOAH compatible

**Options:**

- Carrying case

**AWK-MT10**

Kamlex MT10

**£2995.00**

**AWK-ACC20**

Hard carry case for above

**£216.00**

**Amplivox 116 Screening Audiometer**

- Audiometer of choice for primary care and school use
- Portable and simple to use
- Designed and manufactured in the UK
- Reduced test time with threshold retention function
- 10 comprehensive test frequencies
- -10 to 100 dBHL intensity range
- Supplied with internationally approved headphones (TDH39)
- Includes patient response button, audiometric headset, carry case and concise operating instructions



**AU116** Amplivox 116 Screening Audiometer

**£685.00**

**AU116.3** Battery power function

**£98.00**

**Amplivox Calibration Services:**

All audiometers should be calibrated at regular intervals.

To arrange the calibration of an Amplivox audiometer or tympanometer please call Amplivox on 01865 842411

Your local medical supplier is just a phone call away. See cover for full details  
Prices shown are exclusive of VAT and delivery

**Interacoustics AS608 Audiometer****Ideal for clinics and school nurses**

Portable robust lightweight audiometer suitable for measuring air conduction thresholds in schools, industry and wherever quick hearing evaluations are performed.

- Robust construction
- 11 frequencies 125Hz to 8kHz
- Customise frequency test selections
- Output levels from -10dB to 100dBHL output
- Pure tone, warble tone and pulsed tone
- Battery status indicator
- Carrying bag
- Uses standard AA batteries

**Options:**

- Patient response switch
- Noise reducing headset
- Ext. power supply

**AWK-AS608** Kamlex AS608 Audiometer

**£815.00**

**AWK-215** Patient response switch

**£36.00**

**TS608** Mains adaptor

**£96.00**



Your local medical supplier is just a phone call away. See cover for full details  
Prices shown are exclusive of VAT and delivery

Section B **Audiometry Audiometers****Interacoustics AS216 Audiometer****Ideal for GPs and private hospitals**

Portable lightweight screening audiometer ensuring fast and reliable testing. The manual or automatic test facility on the AS216 when used with the talk forward function makes the AS216 easy to work with, particularly in sound cabin installations.

- Fully portable
- 11 frequencies 125Hz to 8kHz
- -10dB to 120dBHL output
- Standard TDH39 headset
- Bone conductor
- Silent stimulus switch
- Automatic threshold switch
- Pure, warble and pulsed tones
- Built in Talk Forward
- LCD display
- PC on line monitoring available

**Options:**

- Insert phones
- Noise reducing headset
- Carrying case



<b>AWK-AS216</b>	Kamplex AS216 Audiometer	<b>£1465.00</b>
<b>AWK-E5A</b>	E5A insert headphones	<b>£275.00</b>
<b>AWK-ACC26</b>	Hard carry case	<b>£196.00</b>
<b>AWK-KSB25</b>	Soft carry case	<b>£60.00</b>

**Interacoustics AD226 Diagnostic Audiometer****Full featured air/bone audiometer**

Suitable for any stationary or portable application where diagnostic testing of air and bone conduction thresholds is performed.

- 125Hz - 8000Hz
- -10dB to 120dBHL output
- Air and bone conduction
- WN and NB masking
- 5dB and 1dB attenuators
- Pure, warble and pulsed tones
- Automatic threshold test
- SISI and Bekesy tests
- Stenger and ABLB tests
- Tone in noise test (Langenbeck)
- Talk forward
- USB/RS232 interfaces
- PC on line monitoring available
- NOAH compatible
- Printing options via PC

**Options:**

- Insert phones
- Noise reducing headset
- Carrying case



<b>AWK-AD226</b>	<b>£1695.00</b>	<b>AWK-ACC26</b>	<b>£196.00</b>
Kamplex AD226 audiometer		Hard carry case	
<b>AWK-E5A</b>	<b>£275.00</b>	<b>AWK-KSB25</b>	<b>£60.00</b>
E5A insert headphones		Soft carry case	

Your local medical supplier is just a phone call away. See cover for full details  
Prices shown are exclusive of VAT and delivery

Headsets / Tuning Forks **Audiometry Section B****Audiocups and Headsets**

- Noise reducing enclosures to fit over headsets and reduce ambient noise when testing

Image shows enclosures fitted over headset



**AU601** Audiocups - noise reducing enclosures **£98.00**

**THD39** Replacement headsets **£210.00**

Please note that replacement headsets must be calibrated to each machine and it is therefore necessary for your audiometer to be returned for calibration and fitting of headset.

**AW Spirit Tuning Forks**

- Lightweight alloy tuning forks
- Precision frequency balanced
- Very reasonably priced
- Choice of 128Hz or 256Hz



**G.650.00** AW Spirit 128hz Tuning Fork **£22.00**

**G.650.01** AW Spirit 256hz Tuning Fork **£22.00**

**Tuning Forks****Tuning Forks**

**G.600.00** Gardiner Brown Tuning Fork with foot, 128Hz **£22.95**

**G.600.01** Gardiner Brown Tuning Fork with foot, 256Hz **£20.95**

**G.600.02** Gardiner Brown Tuning Fork with foot, 512Hz **£20.95**

**G.600.03** Gardiner Brown Tuning Fork with foot, 1024Hz **£20.95**

**G.600.04** Gardiner Brown Tuning Fork with foot, 2048Hz **£21.95**

**G.600.05** Gardiner Brown Tuning Fork with foot, 4096Hz **£21.95**

**G.600.10** Gardiner Brown Tuning Fork, set of 6 in box **£172.95**



**36.224.10** Martin Rydel Seiffer Tuning Fork with adjustable frequency **£136.50**

Your local medical supplier is just a phone call away. See cover for full details  
Prices shown are exclusive of VAT and delivery