

Section C Blood & Urine Analysis

Bayer Reagent strips for Urinalysis



Please note:- This only represents a selection of the most popular Bayer reagent strips and tests available. Should the required test not be listed, we will be happy to source it for you. Please call our sales line for further details.

2857	Uristix - 50 test strips Tests performed:- Glucose; Protein	£14.50
2743	Labstix SG - 100 test strips Tests performed:- Glucose; Ketones; Specific gravity; Blood; pH; Protein	£42.00
2741	Multistix SG - 100 test strips Tests performed:- Glucose; Bilirubin; Ketones; Specific gravity; Blood; pH; Protein; Urobilinogen	£48.00
2300	Multistix 10 SG - 100 test strips Tests performed:- Glucose; Bilirubin; Ketones; Specific gravity; Blood; pH; Protein; Urobilinogen; Nitrite; Leucocytes	£56.00
2304	Multistix 8 SG - 100 test strips Tests performed:- Glucose; Ketones; Specific gravity; Blood; pH; Protein; Nitrite; Leucocytes	£48.00
2283	Multistix GP - 25 test strips Tests performed:- Glucose; Ketones; Specific gravity; Blood; pH; Protein; Nitrite; Leucocytes	£14.50
2872	Albustix - 50 test strips Tests performed:- Protein	£10.50

CLINITEK status Urine Chemistry Analyser

- The efficient way to get dependable results in any patient setting
- User friendly with display prompts, making the Clinitek status intuitive and easy to use
- Requires only 10 seconds of operators time, full patient results printed in 1 minute
- Abnormal results are flagged for quick review and consideration of treatment
- Clinitek status tests the following:- Leucocytes; Nitrite; Urobilinogen; Protein; pH; Blood; Specific Gravity; Ketone; Bilirubin; Glucose
- Screens out non-infected urine samples so that only positives are referred for laboratory follow-up
- 3 step process:
 - 1 Dip reagent strip into sample and press start button
 - 2 Blot reagent strip and place on instrument feed table
 - 3 Wait one minute for results to be printed
- Uses Multistix 10SG, Multistix 8SG or Multistix GP reagent strips
- Automatic self-tests ensure reliable and accurate results
- Audible alarm indicates errors immediately
- Prints up to 3 copies of results for patient records
- 12 month warranty



5234598	CLINITEK status Urine Chemistry Analyser	£795.00
8727986	Thermal printer rolls, pack of 5	£39.00

Your local medical supplier is just a phone call away. See cover for full details
Prices shown are exclusive of VAT and delivery

Blood & Urine Analysis Section C

Bayer Glucolet® 2 with Minilet™ Lancet



- Glucolet 2 and Minilet systems provide fast, safe and simple blood sampling
- Glucolet 2 is easy to use:
 - 1 Press the transparent barrel to click into position;
 - 2 Fit a new Minilet Lancet;
 - 3 Position over the patients fingertip;
 - 4 Press the release button
- Greater safety from cross infection - The Glucolet 2 can only be re-primed once the used Minilet has been removed for disposal
- Minilet lancet is only exposed the instant the patients finger is pricked, protecting against accidental needle-stick injury
- New superfine Minilet Lancet for improved comfort
- Strictly single use only Minilets



5977	Glucolet® 2	£18.00
5971/200	Minilet Lancets, box of 200 (8 x 25)	£49.00
5971/25	Minilet Lancets, box of 25	£7.00

Blood Lancets


- Stainless steel
- Sterile
- Economically priced



1110101	£7.50
Box of 200 single use, sterile, blood lancets	

Safety Blood Lancets

- 4 different needle sizes available
- Colour coded
- Available in 4 sizes
 - I. 26g with 1.8mm depth - blood glucose testing (Yellow)
 - II. 21g with 1.8mm depth - thin skin (Orange)
 - III. 21g with 2.4mm depth - normal skin thickness (Pink)
 - IV. 18g with 1.8mm depth - for baby's heel (Green)
- Economic box of 100



1110115	26G with 1.8mm depth (Yellow)	£19.99
1110116	21G with 1.8mm depth (Orange)	£19.99
1110117	21G with 2.4mm depth (Pink)	£19.99
1110118	18G with 1.8mm depth (Green)	£19.99

Allergy Lancets

- Specially designed small tip for allergy testing
- 0.9mm tip length
- Stainless steel, single packed, sterile
- Economic box of 200



1110106	£11.50
Allergy lancets with 0.9mm tip	

Lancing Device

- For painless blood sampling
- Easy regulation of penetration depth
- Use with 1110104 Blood Lancets



1110110	Lancing device for use with 1110104 lancets	£11.00
1110104	Sterile blood lancets for use with above pack of 100	£4.75

Your local medical supplier is just a phone call away. See cover for full details
Prices shown are exclusive of VAT and delivery