

Section K Home Care Therapy

Hippo Humidifier - Steam Vaporiser

- Large capacity environmental humidifier
- Relieves bronchial disorders and promotes good respiratory health
- Stable and robust build quality
- The 5.5 litre tank is moulded from tough ABS plastic with hand-holds for manoeuvring
- Pull-out tray for aromatherapy oils
- External water level indicator
- Steam flow adjustment
- Safe to use with fused cut-out, complies with EN 60 335-2-98/A1 standard
- Up to 12 hours continuous use
- Substantial element system ensures the longevity of the unit



120.50.000 Hippo Humidifier £54.00

AW Medex Massager

- Large surface area massager for effective relief from muscular strain and injury
- Effective with all large muscle groups, especially back injuries
- Used by professionals and private patients alike
- Mains 240v powered
- Four-point vibration control ensures equal vibration across entire surface of pad
- Includes 2 hardwearing rubber pads, 1 x smooth and 1 x hyperaemisation
- Supplied in smart carrying case



100.10.000 AW Medex large surface area massager £175.00



boso
BOSCH + SOHN

BoSo Relax 8000

- High quality personal massager
- Anatomically shaped handle for easy application
- Rechargeable
- Extendable to reach awkward areas
- Multiple-surface areas to target specific regions

BoSo 8000 BoSo 8000 Massager £59.50

Epicondylitis Clamps

- Safe and effective relief from the pain of tennis elbow
- Lightweight and comfortable to wear
- Allows instant relief and long-term healing
- Available in a large range of sizes



105.20.225	Epicondylitis Clamp 225mm - X-Small	£21.50
105.20.245	Epicondylitis Clamp 245mm - Small	£21.50
105.20.265	Epicondylitis Clamp 265mm - Med	£21.50
105.20.290	Epicondylitis Clamp 290mm - Large	£21.50
105.20.325	Epicondylitis Clamp 325mm - X-Large	£21.50

Therapy Home Care Section K

Light therapy

Seasonal Affective Disorder (SAD) or the "Winter Blues", effects up to 30% of the UK population to varying degrees.

Light therapy has been consistently proven to alleviate seasonal symptoms in approximately 85-90% of cases, within 4-12 days of use.

There are four classic symptoms experienced by SAD sufferers during the winter months:

- Fatigue and lack of energy
- Greater need for sleep/sleep disorders
- Cravings for carbohydrates and sweet foods, often leading to weight gain
- Depression, low moods

The Brite-Box Sun Stimulator 410D-Lux

- Delivers 10,000 Lux (6,000K)
- Certified 93/42/EEC Annex V for medical devices
- Mains operated 4x36W flicker free compact lamps
- Low EMF and no harmful UV emissions
- Lightweight and portable
- 30 minute timer
- Polycarbonate screens and durable ABS casing
- Average 10,000 hours tube life
- Ideal for use within the home or office



Brite-Box Brite-Box Sun Stimulator 410D-Lux £199.00

**BoSo Infrared Therapy Lamps**

boso
BOSCH + SOHN



- 100 or 150 watt IR bulb
- Functional design with stable platform
- 240v mains operation

BoSotherm 4000	100w I.R. Therapy lamp	£19.95
BoSotherm 4100	150w I.R. Therapy lamp	£34.00

BoSo Pocket Laser 9000

boso
BOSCH + SOHN



- Active low output laser ideal for minor skin irritations
- Battery powered

BoSo 9000 BoSo Pocket Laser £78.00

BoSo Acupuncture 9100

boso
BOSCH + SOHN



- Needle free acupuncture
- Includes rechargeable batteries

BoSo 9100 BoSo Acupuncture Laser £165.00

Your local medical supplier is just a phone call away. See cover for full details
Prices shown are exclusive of VAT and delivery

Your local medical supplier is just a phone call away. See cover for full details
Prices shown are exclusive of VAT and delivery

Section K Home Care Pain Relief

Retro Tens

The **DX6605D** is a modern, state of the art TENS designed with input from therapists and their patients, providing symptomatic relief of both chronic and acute pain.

- Dual Channel TENS machine with a distinctive retro feel with its bright chrome controls and chrome bezel ring
- Large clear LCD display provides the modern touch and shows the selected programs, frequencies and times
- High quality electrodes and leads
- Clear instructions for ease of use
- For the beginner new to TENS, 12 pre-set programs provide a wide choice of therapy
- Advanced user can select from 2 user defined programs with the aid of his medical advisor and tailor programs to their individual need
- Attractive and discrete packaging with optional belt pack available



DX6605D		£39.00
Retro Twin Channel Tens with 4 gel pads and connecting leads		
DX66	Gel Pads (pair)	£1.75 per pair
DX65	Electrode cable	£7.00
DX64	Optional Belt pouch	£7.00

Med-Sport Dual Channel Tens



Attractive yet discretely presented TENS therapy device, the DX6605E is a modern, state of the art TENS designed with input from therapists and their patients and it provides symptomatic relief of both chronic and acute pain.

- Soft touch exterior is tactile yet durable
- Chrome controls and clear display provide quality visual appeal
- The DX6605E unit offers a choice of 13 preset programs including sports modes for muscle toning and warm up
- Very simple to use the DX6605E is supplied with high quality electrodes and leads
- Optional belt pouch is available

DX6605E		£43.00
Med-Sport Twin Channel Tens with 4 gel pads and connecting leads		
DX66	Gel Pads (pair)	£1.75 per pair
DX65	Electrode cable	£7.00
DX64	Optional Belt pouch	£7.00

AW Long Handled Nail Scissors

- Serrated nail scissors designed to cut tough nails that are difficult to reach
- Toe nails can be reached easier
- Serrated blades make light work of the thickest nails
- Made from stainless steel with autoclaveable plastic handles
- 20cm overall length
- Suitable for the elderly




99.05.070	AW Long Handled Nail Scissors	£5.95
------------------	-------------------------------	--------------

Your local medical supplier is just a phone call away. See cover for full details
Prices shown are exclusive of VAT and delivery

Pillows, Pads, Warmers Home Care Section K

BoSo Heat Pads

A range of high quality heat pads ideal for localized heating, available in a range of sizes.



BoSoTherm 1000	30x40cm heat pad with 2 heat settings. 50w	£26.00
BoSoTherm 1100	30x40cm heat pad with 3 heat settings. 50w	£34.00
BoSoTherm 1200	30x40cm heat pad with timer control. 50w	£42.00

BoSo Foot Warmer



BoSo Foot Warmer	2 heat position foot warmer. 24w	£42.00
-------------------------	----------------------------------	---------------

BoSo Heated Blankets



BoSo 2000	130 x 75cm, 3 temperature positions. 50w	£38.00
BoSo 2100	150 x 80cm, large size. 60w	£48.00
BoSo 2200	150 x 80cm, large size with 90 minute timer. 60w	£58.00

BoSo Back Pad



BoSo Heated Back Pad	Heated back pad with 60 minute timer. 56w	£44.00
-----------------------------	---	---------------

Our team of advisors will be happy to discuss your patients requirements with yourself or your patients directly. For a professional and confidential Home Care service, call the usual number at the front of the catalogue.


PRO-PIL-O

PRO-PIL-O Vario:

Internationally known neck support pillow. Soft contoured special foam keeps the head, neck and shoulders in a relaxed position during sleep.



PPV1	£39.00
-------------	---------------



PRO-PIL-O Traveller:

Traveller gives continued sleep support, on holiday flying or travelling.

PPT1	£35.00
-------------	---------------

AW Grant Dyna-Care Ripple Mattress

- Straight forward pressure relief, at a reasonable price
- Gives uniform patient pressure relief, for any single bed
- 21 individual pressure cells inflate and deflate on an alternating basis, ensuring pressure equalisation at every point
- Fits under the base sheet of any single bed
- Ripple control unit operates virtually silently, and there is no noise or vibration within the mattress, allowing the patient a sound nights sleep, in comfort
- 2 year warranty



102.00.000	AW Grant Dyna-Care Ripple Mattress, complete	£194.00
-------------------	--	----------------

Your local medical supplier is just a phone call away. See cover for full details
Prices shown are exclusive of VAT and delivery

Section K Home Care Crutches

Kowsky Crutches

A brand new range of child and adult crutches, designed for the long-term user and made to the highest standard, in Germany. Lightweight, aluminium design with many quality features:

- adjustable length closed-cuff forearm band (Vario)
- hinged-cuff allows crutch to hang loose when not in use (Vario)
- adjustable height with external pin and fixing collar
- wide choice of bright colours available
- ergonomic and anatomic handles, available in hard or soft-feel
- colour coded accessories
- 'clatter' free - no loose rattle noise associated with normal crutches

Vario Adult Crutches

Adult Vario crutches are supplied with black ABS plastic arm and handle assemblies, combined with a choice of the following colour finishes:

- Colours available:
 1. Red
 2. Blue
 3. Turquoise
 4. Purple
 5. Black
 6. Yellow
 7. Blackberry
 8. Bronze
 9. Aluminium



variable
and
extendable



Vario – adult:

* Please add the number of the colour required

BWC-KVA01*	£49.90
Vario adult anatomic soft-feel grips, pair	
BWC-KVA02*	£58.50
Vario adult anatomic normal grips, pair	
BWC-KVA03*	£28.50
Vario adult ergonomic soft-feel grips, pair	
BWC-KVA04*	£37.90
Vario adult ergonomic normal grips, pair	

Crutches Home Care Section K

Vario Child Crutches

All the features of the Vario adult, plus:

- 8 individual colour combinations
- Open cuff design for convenience
- Available with ergonomic or anatomic soft-feel grips



Vario – child:

* Please add the number of the colour required

BWC-KVC01*	£49.90
Vario child anatomic soft-feel grips, pair	
BWC-KVC02*	£41.90
Vario child ergonomic soft-feel grips, pair	

Vario Rubber Crutch Tips

Are available in a wide selection of bright colours, with steel reinforcement and made to a very high standard for longevity.



BWC-KF01	£1.40
Crutch tip, grey	
BWC-KF02	£1.40
Crutch tip, black	
BWC-KF03	£1.40
Crutch tip, brown	
BWC-KF04	£1.40
Crutch tip, red	
BWC-KF05	£1.40
Crutch tip, blue	
BWC-KF06	£1.40
Crutch tip, yellow	
BWC-KF07	£1.40
Crutch tip, turquoise	
BWC-KF08	£1.40
Crutch tip, purple	
BWC-KF09	£1.40
Crutch tip, blackberry	

Tips are available in both child and adult sizes.
Please add 'A' to the end of the order code for Adult or 'C' for Child.

Your local medical supplier is just a phone call away. See cover for full details
Prices shown are exclusive of VAT and delivery

Your local medical supplier is just a phone call away. See cover for full details
Prices shown are exclusive of VAT and delivery