


Compact Quality without Compromise

MyLabTM15





No Compromises

- > Premium Black & White Image Quality
- > Upgradeable to Color and Doppler
- > Three Transducer Ports for Convex, Linear and Endocavity Probes

MyLabTM15



> MyLab World Advanced Technology

The MyLab15 is Esaote's new compact, premium black & white ultrasound system, which can be upgraded to include color and Doppler. Engineered with Esaote's ESA-DYN and XView technologies, the MyLab15 image quality is truly exceptional for a system in its size and price class. XView is the next generation of Esaote's proprietary solution for superior contrast resolution. The ESA-DYN microchip allows Esaote to offer the MyLab15's exclusive high-dynamic signal processing.

> Improved Workflow through Design

The MyLab15's ergonomic design defines it as a system created for ease-of-use and increased practice workflow and efficiency. Every system has three probe connectors and two rest ports to be used with various high-density convex and linear transducers, including Esaote's unique 198° field-of-view endocavity probe. This eliminates the time needed to change probes during or between examinations. Additionally, the MyLab15 is available with either a SVGA monitor or LCD screen depending on the users' preference.





> Innovative USB Solutions

The MyLab15 also integrates today's growing USB port functionality. Every system is supplied with a personal USB memory drive, and the built-in USB port allows the user to save images quickly and easily, as well as connect a PC printer directly to the MyLab15. The memory drive can be used in any PC to save, print or further process the images. A perfect solution for today's busy OB/GYN practice: patients can be given a USB memory drive to bring to each visit to take copies of their images home to share.

- > Modern USB Image Storage and Transfer
- > Expand Patient Satisfaction



- > General Imaging and OB/GYN Application Modules
- > Mother-Child Database
- > Upgradeable to 4D

MyLab15



> Complete Application & Feature Possibilities

Based on the MyLab World modular architecture, the MyLab15 offers a wide range of application packages and can be individualized to every user's specific needs. The General Imaging and OB/GYN modules include: Abdomen, Musculoskeletal, Thyroid, Breast, Obstetrics and Small Parts. The MyLab15's modularity also includes special features, such as TEI™ (Tissue Enhancement Imaging). The easy-to-use TEI feature, with its one-button access and quick response time, is a true technological advantage in routine ultrasound use. The MyLab15 is also upgrade-able to Esaote's revolutionary 4D imaging configuration, which truly adds a new dimension to obstetric examinations.

> Added OB/GYN Benefits

The MyLab15 OB/GYN module contains a new feature specifically designed for the OB/GYN ultrasound user. The Mother-Child Database is an internal database system that allows the clinician to save fetal measurements made at each visit and to view and analyze these over time, all in one easy report. The Mother-Child Database feature provides the MyLab15 OB/GYN user with an extra tool for providing the best diagnostic care possible.



Powerful



MyLabTM15



0344

ISO 13485: 1996

ISO 9001: 2000

www.esaote.com

ESAOTE S.p.A.

Headquarters: Via Siffredi, 58 16153 Genova, ITALY, Phone +39-010-6547.1, Fax +39-010-6547.275
International Activities: Via di Caciolle, 15 50127 Firenze, ITALY, Phone +39-055-4229.1, Fax +39-055-4229.208, E-mail: international.sales@esaote.com

**France
ESAOTE
France S.A.R.L.**

22, rue Pierre Grange
94124 Fontenay-sous-Bois
Cedex
Phone +33-1-4871.2525
Fax +33-1-4871.3630
E-mail: esaote.france@wanadoo.fr

**Germany
ESAOTE BIOMEDICA
Deutschland GmbH**

Hanns-Braun-Straße 50
85375 Neufahrn bei Freising
Phone +49 180 5372683
Fax +49 8165 61820
E-mail: esaote@esaote.de

**Russian Federation
and CIS
ESAOTE S.p.A.**

Leningradskiy prospect
bld. 18, office 5 and 6.
125040 Moscow, Russia
Phone +7-095-232.18.33
Fax +7-095-232.02.05
E-mail: esaotemoscow@mtu-net.ru

**Spain
ESAOTE
ESPAÑA S.A.**

Avda San Sebastian, s/n
08960 Sant Just Desvern
Barcelona Spain
Phone +34 93 473 2090
Fax +34 93 473 2042
E-mail: info@esaote.es

**USA
BIO SOUND ESAOTE Inc.**

8000 Castleway Drive
P.O. Box 50858
Indianapolis
Indiana 46250
Phone +1-317-813.6000
Fax +1-317-813.6600
E-mail: biosound@biosound.com

**China
ESAOTE CHINA Ltd**

Unit 1903-5, 19/F COSCO Tower
Grand Millennium Plaza
No.181-183 Queen's Road Central
Hong Kong
Phone +852-2545.8386
Fax +852-2543.3068
E-mail: esaote@esaotechina.com

**The Netherlands
PIE MEDICAL Benelux B.V.**

P.O. Box 1132
6201 BC Maastricht
The Netherlands
Phone +31-43-3824650
Fax +31-43-3824651
E-mail: benelux@pie.nl